

D.B.M. C-Lock kateterlås

D.B.M. C-Lock 30 % & 46,7 % er et CE-merket natriumtricitrat-kateterlåsemiddel.

D.B.M. C-Lock holder sentrale venekatetre (CVC) åpne mellom behandlinger, samtidig som det beskytter kateteret mot mikrobiell vekst og reduserer dannelsen av biofilm.

Produktet er pakket i sterile poser med to sprøyter/pose.



Bestillingsinformasjon

Bestillingsnummer	Produktnavn	Pakning
PCE15KA2	C-Lock 30 % 2 x 2,5 ml	120 st
PCE14KA2	C-Lock 46,7 % 2 x 2,5 ml	120 st

Referens från USA:

1. Weijmer, M.C., van den Dorpel, M.A., Van de Ven, P.J.g., ter Wee, P.M., van Geelen, J.A. C.A., Groeneveld, J.O., van Jaarsveld, B.C., Koopmans, M.G., Poole, C.Y., Schrander Van der Meer, A.M., Siegert, C., & Stas, K.J.F. (2005). Randomized, clinical trial comparison of trisodiumcitrate 30% and heparin as catheter-locking solution in hemodialysis patients. Journal of American Society of Nephrology, 16, 2769-2777.
2. Ash, S.R., Mankus, R.A., Sutton, J.M., Criswell, R.E., Crull, C.C., Velasquez, K.A., Smeltzer, B.D., & Ing, T.S. (2000). Concentrated sodium citrate (23%) for catheter lock. Hemodialysis International, 4, 22-31.
3. Winnett, G., Nolan, J., Miller, M., & Ashman, N. (2008). Trisodium citrate 46,7% selectively and safely reduces staphylococcal catheter-related bacteremia. Nephrology Dialysis Transplant, 23,3592-3598.
4. Lok, Charmin E., Debra Appleton, Cynthia Bhola, Brian Khoo, and Robert M.A. Richardson "Trisodium citrate 4% - an alternative to heparincapping of haemodialysis catheters." Nephrology Dialysis Transplantation. 22. (2007): 477-483. Print.

Spesifikasjoner

- Sterile 2,5 ml ferdigfylte sprøyter med luer-hette-tilkobling gir enkel tilkobling og støtter aseptisk arbeid
- Gjennomsiktige polypropylenssprøyter av høy kvalitet, med en lateksfri stempelgummitetning.
- Produktinformasjon og måleskala tydelig på siden av sprøyten. Fargekodede korker for ulike styrker
- Laget i Italia i henhold til ISO 13485 for kvalitet og pålitelighet.
- Produktets emballasje muliggjør god holdbarhet og enkel oppbevaring ved en temperatur på +5 °C til +30 °C.
- Gyldighetstid 3 år i uåpnet emballasje.