

ClearGuard™ HD

antimikrobiellt skyddande hölje för hemodialyskatetrar

ClearGuard HD antimikrobiell skyddskork för hemodialyskateter skyddar även hemodialyskateterns lumen mellan behandlingarna.

ClearGuard HD skyddar kateterlumen både externt och invärtes med hjälp av sin klorhexidinstav, vilket har visat sig minska risken för kateterrelaterade infektioner

- ClearGuard HD antimikrobiell hemodialyskateterlock sitter på mellan behandlingar.
- Kateterförslutningen håller det antimikrobiella medlet inne i katetern i lumen mellan behandlingarna.
- ClearGuard HD-kåpor används istället för ett vanligt skyddskork



- Latex fri
- Gamma steriliserad
- 2 st/förpackning
- Max drifttid 72 timmar

Beställningsinformation

Produktkod	Produkt	Förpackning
011-CGHD-100	ClearGuard	100 st



National
Kidney
Foundation®

KDOQI
Kidney Disease
Outcomes Quality Initiative

ClearGuard HD caps recommended in NKF's KDOQI Clinical Practice Guideline for Vascular Access: 20195
21.3 KDOQI considers it reasonable to use an antimicrobial barrier cap to help reduce CRBSI in high-risk patients or facilities; the choice of connector should be based on clinician's discretion and best clinical judgment. (Expert Opinion)