

ClearGuard™ HD

antimikrobinen suojakorkki hemodialyysikatetreille

ClearGuard HD -antimikrobinen hemodialyysikatetrin suojakorkki on suunniteltu suojaamaan hemodialyysikatetrin lumenia hoitojen välillä. Ainutlaatuinen antimikrobinen klorhexidin-sauva korkissa suojaa katetrin lumenia sekä ulkoa- että sisältäpäin, mikä tutkimusti pienentää katetriperäisten infektioiden riskiä. Soveltuu käytettäväksi tavanomaisempien hemodialyysikatetrien kanssa.

- ClearGuard HD -antimikrobinen hemodialyysikatetrin suojakorkki työnnetään lukkoaineella täytettyyn katetriin
- Katetrin sulkija pitää antimikrobisen aineen katetrin lumenen sisällä hoitojen välillä
- ClearGuard HD -korkkeja käytetään tavallisen korkin sijasta



- Lateksiton
- Gammasteriloitu
- 2 kpl/pkt
- Maksimi käyttöaika 72 h

Tilaustiedot

Tilausnro	Tuote	Pakkaus
011-CGHD-100	ClearGuard HD-suojakorkki	100 kpl



National
Kidney
Foundation®

KDOQI®
Kidney Disease
Outcomes Quality Initiative

ClearGuard HD caps recommended in NKF's KDOQI Clinical Practice Guideline for Vascular Access: 20195
21.3 KDOQI considers it reasonable to use an antimicrobial barrier cap to help reduce CRBSI in high-risk patients or facilities; the choice of connector should be based on clinician's discretion and best clinical judgment. (Expert Opinion)