

Comfort-3 Flex

toimenpidetuoli

- The Comfort-3 Flex Therapy Chairs are designed to assist and facilitate the work of physicians and medical staff
- They increase the patient's comfort considerably and therefore help to relieve tension during treatment
- Several optional accessories for optimal choice
- Slim design, less space needed



Comfort-3 Flex

ART.-NR. 723C3F01

Optional accessories	Art.-Nr.
EPN kit incl. multicontact pin with internal connecting cable HIGHLY RECOMMENDED FOR DIALYSIS	713E004
Trendelenburg foot-switches on both sides of underframe	713E005
Rechargeable battery back-up operation in case power failure (24V DC)	713E006
Reading lamp with LED light	713E008
USB charging unit for smart devices (5V, 1A)	713E009
Connecting cable (3 m) for potential equalization (EPN)	713E012
Comfort-T2 Multifunctional table	713M009
Flexible height adjustable armrest for universal medical applications (one pair)	713M010
Integrated gas spring in the arm support for smoother and safer operation (one pair)	713M011
Push bar	713M012
Universal holder for tablets and smartphones	713M013
Handle on both sides of the seat	713M014
Motorised height adjustment	713M017
Central operating brake with 4x diam. 7.5 cm twin wheel castors	713M020
4x diameter 10 cm swivel castors with total lock	713M021
Motorized footrest	713M032
Paper roll holder	713M033
Infusion bottle holder (IV pole) with height adjustment	713M038
Manual quick release of backrest for CPR (Cardiopulmonary Resuscitation)	713M083
Individual logo on the upholstery	713U011
Headrest pillow square	713U012
DuoShape upholstery in individual colour combination (40 diff. colours) TUNDRA	713U026
Transparent protecting plastic cover for the footrest	713U027
Transparent protecting plastic cover for the legrest adjustable between 0° and -30° below horizontal	713U028
DuoShape upholstery with 110 mm thickness in standard Chrom-Perle colour combination	713U031
DuoShape upholstery with 110 mm thickness in individual colour combination (40 diff. colours) TUNDRA	713U032
DuoShape upholstery with memory foam in standard Chrom-Perle colour combination	713U033
DuoShape upholstery with memory foam in individual colour combination (40 diff. colours) TUNDRA	713U034
DuoShape upholstery with flame retardant artificial leather combination (15 diff. colours) PALMA	713U035

General technical data			
Fixed seat height	62 cm	Dimensions in seat position	L 145 x W 87.5 cm
Optional height adjustment range	55–70 cm	Size of castors	4 x Ø 7.5 cm twin wheel castors with total lock
Backrest adjusting range	80° (-14° under horizontal in Trendelenburg position)	Weight of Chair	65 kg
Seat tilt adjustment	14°	Maximal load	200 kg weight capacity
Legrest adjustment	30°	Standard Power Supply	100-240 V AC, 50/60 Hz
Footrest adjustment	25 cm	Technical data of motors	24V DC, IPX6
Size of upholstery	L 210 x W 60 x H 9.5 cm	Frame colour	RAL 9001
Armrest size	L 57 x W 15 x H 6.5 cm	Applied directive and standards	(EU) 2017/745 medical device regulation and harmonized standards