

NOW AVAILABLE **X-Series**

PD Cath ORDERING INFORMATION

"X" Series		Description
Catheter Only 5/Box	Sets 5/Box	
MC20XC	MC20XS	15F x 64CM X-Series

"V" Series		Description
Catheter Only 5/Box	Sets 5/Box	
MC200301	MC200101	Straight, 43cm left V-NECK PD Catheter w/ Double Cuff
MC200302	MC200102	Straight, 43cm right V-NECK PD Catheter w/ Double Cuff
MC200401	MC200201	Coiled, 63cm left V-NECK PD Catheter w/ Double Cuff
MC200402	MC200202	Coiled, 63cm right V-NECK PD Catheter w/ Double Cuff

"I" Series		Description
Catheter Only 5/Box	Sets 5/Box	
MDC31S	---	Straight, 31cm PD Catheter w/ Double Cuff
MDC32S	---	Straight, 32.5cm PD Catheter w/ Double Cuff
MDC37S	---	Straight, 37cm PD Catheter w/ Double Cuff
MDC42S	MPD-242	Straight, 42cm PD Catheter w/ Double Cuff
MDC47S	---	Straight, 47cm PD Catheter w/ Double Cuff
MDC57C	MPD-257	Coiled, 57.5 PD Catheter w/ Double Cuff
MDC63C	MPD-263	Coiled, 63cm PD Catheter w/ Double Cuff
MSC31S	---	Straight, 31cm PD Catheter w/ Subcutaneous Single Cuff
MSC37S	---	Straight, 37cm PD Catheter w/ Subcutaneous Single Cuff
MSC42S	---	Straight, 42cm PD Catheter w/ Subcutaneous Single Cuff
MSC46S	MPD-146S	Straight, 46cm PD Catheter w/ Subcutaneous Single Cuff
MPP41S	MPD-141	Straight, 41cm PD Catheter w/ Preperitoneal Single Cuff
MPP46S	MPD-146	Straight, 46cm PD Catheter w/ Preperitoneal Single Cuff
MPP40C	---	Coiled, 40.25cm PD Catheter w/ Preperitoneal Single Cuff
MPP57C	---	Coiled, 57.5cm PD Catheter w/ Preperitoneal Single Cuff
MPP60C	MPD-160	Coiled, 60cm PD Catheter w/ Preperitoneal Single Cuff

	Accessories	Qty
MP217	17F Vascu-Sheath Introducer Set	5/box
MC50	Male Luer Lock End Cap	100/box
MBCA	Barbed luer Lock Catheter Adaptor	20/box
MC51	Beta Clamp	100/box
MCS	Catheter Stylet	5/box
MTS	Tunneling Stylet	10/box
MTC	Titanium Connector	5/box
TRX	X-Series Tunneling Rod	5/box
FSX	X-Series Faller Stylet	5/box

Catheter Only Components:

- (1) Catheter
- (1) Female Luer
- (1) Clamp
- (1) Male Luer End Cap

Component Set:

- (1) Catheter
- (1) Female Luer
- (1) Clamp
- (1) Male Luer End Cap
- (1) Introducer Dilator & Sheath
- (1) J/Flex Guidewire
- (1) Scalpel
- (1) Introducer Needle
- (1) Tunneling Tool
- (6) Gauze Pads
- (1) 10cc Luer Lock Syringe

medCOMP®
PD Cath



Introducing
X-Series



PD Cath



Medical Components, Inc.
1499 Delp Dr., Harleysville, PA 19438, U.S.A.
P: 215.256.4201 F: 215.256.1787
www.medcompnet.com

CE2797
ISO13485
Certified Quality System

EC REP

MPS Medical Product Service GmbH
Borngasse 20
35619 Braunfels
Germany

PN2520 Rev. F 8/19

medcompnet.com/pdcath



PD Cath Inspired By Nature For Better Living

Better living is possible through peritoneal dialysis. Which is why Medcomp offers a complete line of peritoneal dialysis catheters available in multiple designs, lengths, and configurations. The benefits of PD are numerous to patients and providers alike.

Together, with a complete line of insertion and catheter accessories, Medcomp's PD Catheter line provides the solution you and your patients need.

With 78% of your patients eligible, it's time to consider PD.

➤ **Increase Patient Satisfaction**

Patients who choose PD report higher quality of life. These patients benefit from flexible scheduling, fewer needle sticks, less diet restrictions, and an overall reduced burden of kidney disease.

➤ **Lower Total Costs**

In comparison to HD, PD patients cost less. In fact, they cost an average of \$19,000 per-patient-per-year less.

One piece female luer adaptor.

Single or double "fluffy" cuffs allow secure tissue in-growth to prevent infection.

Large strategically placed side holes for increased flow and reduction in incidence of omental occlusion.

Curled end for a more gentle delivery into the body.

Trim Line

medcompnet.com/pdcath

PD Cath

