



I-SERIES/ V-SERIES

Indications for Use: The peritoneal catheter is indicated for acute and chronic peritoneal dialysis.

Contraindications: Infected anterior abdominal wall. Patients with extensive intra-abdominal adhesions. Conditions resulting in diaphragmatic tears or abdominal-chest communication. Unresolved peritonitis. Severe respiratory insufficiency: distending the abdomen with large amounts of dialysate may further compromise pulmonary function. The catheter is not intended for the Tenckhoff trocar method of insertion, nor is it meant to be used for any purpose other than indicated. Recent intestinal surgery.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.



Silicone
Material


**Straight,
Coiled**
Tip Designs

15F
French Sizes

**Single,
Double**
Cuffs



COILED TIP



SINGLE OR DOUBLE CUFF

ONE-PIECE FEMALE
LUER ADAPTER**CATHETER ONLY - STRAIGHT**

MDC31S	31CM I-SERIES, STRAIGHT WITH DOUBLE CUFF	5/BOX
MDC32S	32.25CM I-SERIES, STRAIGHT WITH DOUBLE CUFF	5/BOX
MDC37S	37CM I-SERIES, STRAIGHT WITH DOUBLE CUFF	5/BOX
MDC42S	42CM I-SERIES, STRAIGHT WITH DOUBLE CUFF	5/BOX
MDC47S	47CM I-SERIES, STRAIGHT WITH DOUBLE CUFF	5/BOX
MSC31S	31CM I-SERIES, STRAIGHT WITH SUBCUTANEOUS SINGLE CUFF	5/BOX
MSC37S	37CM I-SERIES, STRAIGHT WITH SUBCUTANEOUS SINGLE CUFF	5/BOX
MSC42S	42CM I-SERIES, STRAIGHT WITH SUBCUTANEOUS SINGLE CUFF	5/BOX
MSC46S	46CM I-SERIES, STRAIGHT WITH SUBCUTANEOUS SINGLE CUFF	5/BOX
MPP41S	41CM I-SERIES, STRAIGHT WITH PREPERITONEAL SINGLE CUFF	5/BOX
MPP46S	46CM I-SERIES, STRAIGHT WITH PREPERITONEAL SINGLE CUFF	5/BOX

CATHETER ONLY CONTAINS:

- (1) Catheter
- (1) Female Luer
- (1) Clamp
- (1) End Cap

BASIC SET CONTAINS:

- (1) Catheter
- (1) Female Luer
- (1) Clamp
- (1) Introducer Needle
- (1) Guidewire
- (1) Scalpel
- (1) Peelable Sheath Introducer
- (1) Tunneler
- (1) 10cc Syringe*
- (1) End Cap
- (6) 4" x 4" Gauze

* - Except MPD-118

BASIC SET - STRAIGHT

MPD-237	37CM I-SERIES, STRAIGHT WITH DOUBLE CUFF	5/BOX
MPD-242	42CM I-SERIES, STRAIGHT WITH DOUBLE CUFF	5/BOX
MPD-146S	46CM I-SERIES, STRAIGHT WITH SUBCUTANEOUS SINGLE CUFF	5/BOX
MPD-118	18CM I-SERIES, STRAIGHT WITH PREPERITONEAL SINGLE CUFF	5/BOX
MPD-141	41CM I-SERIES, STRAIGHT WITH PREPERITONEAL SINGLE CUFF	5/BOX
MPD-146	46CM I-SERIES, STRAIGHT WITH PREPERITONEAL SINGLE CUFF	5/BOX

CATHETER ONLY - COILED

MDC57C	57.5CM I-SERIES, COILED WITH DOUBLE CUFF	5/BOX
MDC63C	63CM I-SERIES, COILED WITH DOUBLE CUFF	5/BOX
MPP40C	40.25CM I-SERIES, COILED WITH PREPERITONEAL SINGLE CUFF	5/BOX
MPP57C	57.5CM I-SERIES, COILED WITH PREPERITONEAL SINGLE CUFF	5/BOX
MPP60C	60CM I-SERIES, COILED WITH PREPERITONEAL SINGLE CUFF	5/BOX

BASIC SET - COILED

MPD-257	57.5CM I-SERIES, COILED WITH DOUBLE CUFF	5/BOX
MPD-263	63CM I-SERIES, COILED WITH DOUBLE CUFF	5/BOX
MPD-160	60CM I-SERIES, COILED WITH PREPERITONEAL SINGLE CUFF	5/BOX

ACCESSORIES

MP217	17F VASCU-SHEATH® INTRODUCER SET	5/BOX
MC50	MALE LUER LOCK END CAP	100/BOX
MBCA	BARBED LUER LOCK CATHETER ADAPTOR	20/BOX
MCS	CATHETER STYLET	5/BOX
MTS	TUNNELING STYLET	10/BOX
MTC	TITANIUM BARB CONNECTOR	5/BOX

Refer to the Table of Contents for Important Risk Information regarding this device.